

Ürün/Primer ambalaj	Ürün kodu	Lot numarası	Son kullanım tarihi
Yeast Extract Agar / 60 mm Ø Petri	YEA-60-P	130423-09	13/07/2023

Fiziksel kalite kontrol		
Parametre	Spesifikasyon	Sonuç
Görünüm	Saydam	Uygun
Renk	Sarımsı kahverengi	Uygun
pH (25°C)	6.3 - 6.7	Uygun

Mikrobiyal kontaminasyon		
İnkübasyon şartları	Spesifikasyon	Sonuç
36 ± 2 °C; 44 ± 4 saat	Mikrobiyal çoğalma gözlenmemeli	Uygun

Verimlilik				
Referans besiyeri	-			
Kontrol yöntemi	Kantitatif			
Kontrol suşu	İnkübasyon şartları	Spesifikasyon	Karakteristik reaksiyon(lar)	Sonuç
<i>Escherichia coli</i> ATCC 25922	36 ± 2 °C; 44 ± 4 saat	PR ≥ 0,7	-	Uygun
<i>Bacillus subtilis subsp. Spizizenii</i> ATCC 6633	36 ± 2 °C; 44 ± 4 saat	PR ≥ 0,7	-	Uygun

Lot kalite kontrol onayı	
Tarih	Onaylayan
17/04/2023	Şibanur KOÇ Kalite Kontrol Sorumlusu

#### Referanslar

ISO 11133:2014, Microbiology of food, animal feed and water - Preparation, production, storage and performance testing of culture media

CLSI M22, Quality Control for Commercially Prepared Microbiological Culture Media, 3rd Edition

Australian Society for Microbiology, Guidelines for assuring quality of food and water microbiological culture media. Second Edition, 2014

