

Helico MED

H.pylori Rapid Urease Test

Helicobacter pylori is an important infectious agent that causes stomach ulcers. Therapies that eliminate this microorganism, treat successfully the ulcers. *Helicobacter pylori* produces a very strong and rapid urease enzyme. The identification of urease activity in gastrointestinal mucosa very strongly indicates the presence of this ulcer causing agent, in this tissue.

HelicoMed is a test medium that identifies urease activity in stomach biopsy specimens. Thanks to its new improved original formula it works faster than all other urease tests. **HelicoMed** provides the optimum conditions in which urease of *Helicobacter pylori* is most active. **HelicoMed** is as specific and sensitive as molecular diagnostic methods.

FEATURE

- Rapid Result
- Easy to Use
- Durable
- Storage at room temperature
- Accurate
- Low Cost
- Dry Test advantage



ADVANTAGE

- Positive and negative results in 2 hours
- An extremely simple use
- Durable for 36 months
- No refrigerator required
- Positive control is provided to confirm Negative Results
- Costs much less than other biopsy tests
- No incubation or reagent required



Helico MED⁺

H.pylori Rapid Urease Test

HelicoMed Plus is simply the most effective and easy to use Rapid Urease Test on the market today. An obvious evolution from the old agar gel based tests; **HelicoMed Plus** offers significant improvements that make the task of detecting H.pylori infections in gastric mucosal biopsies faster, easier, less expensive, and more accurate than ever before.

H. pylori is a major cause of diseases of the upper gastrointestinal tract. Eradication of the infection in individuals will improve symptoms including dyspepsia, gastritis and peptic ulcers, and may prevent gastric cancer. While there are several methods of detecting H.pylori, biopsy testing using rapid urease is considered by professionals worldwide and proven in numerous clinical studies to be the "Gold Standard" for detecting active H.pylori infections.

FEATURE

- Rapid Result
- Easy to Use
- Durable
- Storage at room temperature
- Accurate
- Low Cost
- Dry Test advantage



ADVANTAGE

- Positive result in as little as 1 minute, negative result in 1 hour
- An extremely simple use
- Durable for 36 months
- No refrigerator required
- Positive control is provided to confirm Negative Results
- Costs much less than other biopsy tests
- No incubation or reagent required

